

# Speaking Valve




Besmed Speaking Valve helps patients with tracheostomy speak easily. One way valve enables the patient to inhale normally through the tracheostomy. While the patient is exhaling, one way valve will close and air flow pass through the vocal cord to allow the patient to speak.






## Features

- Besmed Speaking valve enables the patients with tracheostomy to be sensitive in the upper airway.
- It helps the patient recover his/her smell, taste and speaking instinct ability.
- This device will accelerate decannulation or expedite ventilator weaning.
- A handy storage box is provided.
- Single patient use.



- a** Speaking Valve (SV-89480)
- b** Tracheostomy Tube
- c** Balloon Deflated

## Order Information

Code	SV-89480	SV-89470	SV-89460
Description	 <p>Speaking Valve for Tracheostomy Tube ID: 15mm</p>	 <p>Speaking Valve with Oxygen Supply Port ID: 15mm</p>	 <p>Speaking Valve for Ventilator Circuit ID: 15mm, OD: 22mm</p>
Package	20pcs/case (including Storage Box)		

INTERACT WITH US



### BESMED HEALTH BUSINESS CORP.

No. 2, Lane 106, Wu-Kong 3rd Rd, Wu-Ku District,  
New Taipei City, Taiwan. 24889

Tel: +886-2-2290-3959 E-mail: [info@besmed.com](mailto:info@besmed.com)  
Fax: +886-2-2299-9076 Web-site: [www.besmed.com](http://www.besmed.com)

